# VEGA TECHNOLOGIES INC. INSTRUCTION MANUAL Medical Nebulizer Kits

Model: A, A01, A02, A03, A04, A05, A06, A07, A08, A09, A10, A11, A12, A13, A14, A15, A17, A18, A19

File Number: MDR-MNK-SPC001

Version/Revision No: 1.0 Date: November 01, 2022



### (Read this instruction manual carefully before use)

### Contents

1.	Introduction	3
2.	Safety precaution	9
3.	Information for operation	.10
4.	Cleaning and maintenance	.11
5.	Service information	.12
6.	Disposal	.12
7. 8	pecifications table	.13
8. 5	ymbols and description	.15
9. C	Contact information	.16

## 1. Introduction

Thank you for purchasing the Medical Nebulizer Kits. It is intended to cooperate with compressor nebulizer to aerosolize physician-prescribed solutions for inhalation.

We encourage you to thoroughly read this guidebook to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

Intended use	The Medical Nebulizer Kits are intended to cooperate with compressor nebulizer to aerosolize physician-prescribed solutions for inhalation.	
Intended User	Professional and lay person (under the guidance of professional)	
Intended patients	Adult, children, infant with respiratory disorder such as asthma, Bronchiectasis, Bronchiolitis, COPD, Cystic Fibrosis, tracheostomy, other Obstructive or Interstitial lung diseases, Post lung transplant, Pulmonary Hypertension, Respiratory Infection.	
CysticFibrosis,chronicobstructivelungIndicationspneumonia,respiratorytractdisease,respiratoryinfection.		
1. This product should not be used on patients, who unconscious or are not breathing spontaneously.Contraindications2. Pregnant or breastfeeding women.3. Patient is allergic to aerosolized drugs.4. Patient with asthma.		

### **1.1 General description**

	<ul><li>5. Pentamidine powder and suspension.</li><li>6. Anesthetics agent.</li></ul>		
Clinical benefit Convey drug droplet for inhalation.			
Intended Environment	This device is intended for use in a medical facility, such as hospital, clinic and doctor office, and in a room of general household.		
Precautions for use	Warnings and cautions described in the instruction manual should be observed		
Solutions used for nebulization	<ul> <li>β2 agonists: Salbutamol/Albuterol Sulfate;</li> <li>Antimuscarinics: Ipratropium Bromide;</li> <li>Corticosteroids: Budesonide, Fluticasone;</li> <li>Combination drug: Ipratropium Bromide and Albuterol</li> <li>Sulfate;</li> <li>Mucoactive Drugs: Normal saline, Hypertonic saline,</li> <li>Dornase alpha /Dnase;</li> </ul>		
	Antibiotics: Colistimethate Sodium.		

### **1.2 Product identification**

		Model
Nebulizer Bottle	i P	CNK-M03
Air Tube	0	CNK-S01

Mouthpiece	all a	CNK-M06
Nosepiece	3500	CUK-M04
Adult Mask	C	CUK-M01
Child Mask	Ì	CUK-M02
Infant Mask	$\bigcirc$	CNK-M14
Connector		CUK-S01

Model	Consist of
	Includes: Nebulizer bottle ( model :
	CNK-M03),
	Air tube (model: CNK-S01),
	Nose piece (model: CUK-M04),
А	Mouth piece (model: CNK-M06),
	Child mask (model: CUK-M02),
	Adult mask (model: CUK-M01) ,
	Infant mask (model: CNK-M14) ,
	Connector (model: CUK-S01)
	Includes : Nebulizer bottle ( model :
	CNK-M03),
	Air tube (model: CNK-S01),
A01	Mouth piece (model: CNK-M06),
	Nose piece (model: CUK-M04),
	Adult mask (model: CUK-M01),
	Child mask (model: CUK-M02)
	Includes : Nebulizer bottle ( model :
	CNK-M03),
A02	Air tube (model: CNK-S01),
	Mouth piece (model: CNK-M06),
	Nose piece (model: CUK-M04),

	Adult mask (model: CUK-M01),		
	Infant mask (model: CNK-M14)		
	Includes : Nebulizer bottle ( model :		
	CNK-M03),		
	Air tube (model: CNK-S01),		
A03			
	Mouth piece (model: CNK-M06),		
	Adult mask (model: CUK-M01), Child mask (model: CUK-M02)		
	Includes : Nebulizer bottle ( model :		
104	CNK-M03),		
A04	Air tube (model: CNK-S01),		
	Mouth piece (model: CNK-M06),		
	Adult mask (model: CUK-M01)		
	Includes : Nebulizer bottle ( model :		
405	CNK-M03),		
A05	Air tube (model: CNK-S01),		
	Mouth piece (model: CNK-M06),		
	Child mask (model: CUK-M02)		
	Includes : Nebulizer bottle ( model :		
A06	CNK-M03),		
	Air tube (model: CNK-S01),		
	Adult mask (model: CUK-M01)		
	Includes : Nebulizer bottle ( model :		
A07	CNK-M03),		
	Air tube(model: CNK-S01),		
	Child mask (model: CUK-M02)		
	Includes : Nebulizer bottle ( model :		
A08	CNK-M03),		
	Mouth piece (model: CNK-M06),		
	Nose piece (model: CUK-M04)		
	Includes : Nebulizer bottle ( model :		
A09	CNK-M03),		
	Air tube (model: CNK-S01),		
	Mouth piece (model: CNK-M06)		
	Includes : Nebulizer bottle ( model :		
A10	CNK-M03),		
	Air tube (model: CNK-S01)		
A11	Includes : Adult mask (model: CUK-M01),		
	Child mask (model: CUK-M02)		
A12	Includes : Nebulizer bottle ( model :		
	CNK-M03)		
A13	Includes : Child mask (model: CUK-M02)		
A14	Includes: Adult mask (model: CUK-M01)		
A15	Includes: Nose piece (model: CUK-M04)		

A17	Includes: Air tube (model: CNK-S01)
A18	Includes: Mouth piece (model: CNK-M06)
A19	Includes: Infant mask (model: CNK-M14)

### 1.3 Description of other device intended to be used in

### combination

Model	CN-01W	CN-01WA	CN-01WB	CN-01WC	CN-01WD
photograph	in m				C.C.
Model	CN-01WG	CN-01WH	CN-01WI	CN-01WK	CN-02WL
photograph					
Model	CN-02WP	CN-02WW	CN-02WZ	CN-02WR	CN-02WU
photograph					
Model	CN-02WX	CN-02WM	CN-02WK	CN-02WG	CN-01WJ
photograph		0.0			
Model	CN-02MX	CN-01WE	CN-01MI	CN-02WA	CN-02WC
photograph			and the second s		

Model	CN-02WD	CN-02WJ	CN-02MB	CN-02MD	CN-02MH
photograph	UNITED OF				
Model	CN-02MJ	CN-02MK	CN-02MN	CN-02MO	CN-02MU
photograph			A CONTRACTOR		Access of the second
Model	CN-02MY	CN-12WA	CN-02WB	CN-02MC	CN-02MR
photograph	e la				C C C
Model	CN-12MA	CN-12ME	CN-01MD	CN-01MF	CN-02WV
photograph	the second secon				
Model	CN-02MF	CN-02MI	CN-02MT	CN-02WN	CN-02ME
photograph					
Model	CN-02MG	CN-02MV	CN-02WE	CN-02WF	CN-02WH
photograph		C,			

Model	CN-02WI	CN-02WO	CN-02WS	CN-02WY	CN-01MN
photograph					•

## 2. Safety precaution

Read all the information in the instruction manual and any other literature included in the box before using the device.

### Marning:

Indicates a potentially hazardous situation which if not avoided, could result in serious injury.

#### (Usage)

- 1) For type, dose, and regime of medication, follow the instructions of your doctor or respiratory therapist.
- 2) The long air tube of this product may cause strangulation of the user if use in wrong way.
- 3) Do not use accessories, detachable parts and materials not described in the instructions for use.
- 4) Do not use in an oxygen rich environment.
- 5) Do not interconnect this device with other equipment not described in the instructions for use.
- 6) If you feel anything unusual during use, please stop using the device immediately and consult your doctor.
- 7) Do not use only water in the Nebulizer Bottle for inhaling purposes.
- 8) Keep the device out of the reach of unsupervised infants and children. The device may contain small pieces that can be swallowed.

- 9) Do not use the device while sleeping or sleepy.
- 10) Store the device in a clean location, make sure that the device is clean before use.
- 11) Do not store the Air Tube while there is moisture or medication remaining inside it.
- 12) Do not use the device where it may be exposed to flammable gas.
- 13) Always dispose of any remaining medication after use and only use fresh medication each time.
- 14) Do not leave the device unattended with infants or person who cannot express their consent.
- 15) Do not reuse on more than one patient.

#### **General safety precautions:**

1) Inspect the device and parts before using them each time, and check that there are not problems. Inparticular, be sure to check the following:

-The Nebulizer Bottle and Air Tube are not damaged.

2) Operate the device only as intended. Do not use the device for any other purpose.


## 3. Information for operation

#### **3.1 Installation and operation procedure**

**Important**: Prior to initial operation, the device should be thoroughly cleaned referring to the "Cleaning and Maintenance" in the manual.

- 1) Gently twist the top cover of the Nebulizer Bottle counter-clockwise to disassemble the Nebulizer Bottle.(Fig.1)
- 2) Fill the cup base of the Nebulizer Bottle with the medication prescribed by your physician.(Fig.2)

- 3) Gently twist the top cover clockwise to re-assemble the Nebulizer Bottle. Be sure the two sections fit well.(Fig.3)
- 4) Attach one end of Air Tube to the cup base of the Nebulizer Bottle.(Fig.4)



- 5) Attach the other end of Air Tube to the tube connector located on the front of the Compressor Nebulizer.
- Attach one of inhalation device per your choice to the top cover of Nebulizer Bottle, but Infant Mask need to add Connector between mask and Nebulizer Bottle.

Remark: Inhalation device is Mouthpiece, Nosepiece, Adult Mask, Child Mask, Infant Mask.

7) Power on Compressor Nebulizer to begin your prescribed treatment.

## 4. Cleaning and maintenance

#### **Cleaning:**

It is recommended to clean the Medical Nebulizer Kits after each use unless there are other instructions by your doctor.

- 1) Rinse the device with pure water;
- 2) Wipe dry the device with clean towel;
- 3) Ensure the device without residual water, then store it in a clean location.

## 5. Service information

Any abuse, operation other than the intended use of this product as outlined in the manual, negligence, accident, or repairs by someone other than a VEGA' s Authorized Service Professional, shall immediately void this warranty.

Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### 6. Disposal

The device and combination devices may come into contact with infectious material and be contaminated during their lifetime. For this reason, the device and its combination devices should be decontaminated before disposal or transportation.

Dispose of the device properly at the end of its service life. According to the European Regulation (EU) 2017/745 the device should not be disposed of with unsorted domestic waste. Look for a careful separation of materials. Consider the country-specific laws and regulations that apply to the disposal of the device. The proper disposal of the device prevents environmental and human damage.

Outside the EU, the appropriate national regulations for disposal must be observed.

### 7. Specifications table

Device name	Medical Nebulizer Kits	
	Nosepiece: CUK-M04	
	Mouthpiece: CNK-M06	
Consist of	Child Mask: CUK-M02	
	Adult Mask: CUK-M01	
	Nebulizer Bottle: CNK-M03	

		Air Tube: CNK-S01		
		Infant Mask: CNK-M14		
		Connector: CUK-S01		
Type of nebulizer		Compressor nebulizer		
Physical property		Refer to MSDS report of raw material in Attachment 3 of folder 01 Device description and specification including variants and accessories		
Chemical property		Refer to MSDS report of raw material in Attachment 3 of folder 01 Device description and specification including variants and accessories		
Performance		MMAD: 5 um		
Recommended fill volume		2 to 5 ml		
Nebulization rate		≥0. 20 ml/min. (0.9% salt water solution)		
Applied part		Type BF Mouthpiece Nosepiece Adult/child/ Infant mask		
Device Lifetime		300 times use		
Device shelf-life		One year		
	Ambient temperature	5℃~40℃ (41°F~104°F)		
Operating environment	Relative humidity	15%~90% 70КРа~106КРа		
	Atmospheric pressure			
	Ambient temperature	-25℃~70℃ (-13℉~158℉)		
Storage environment	Relative humidity	0%~90%		
	Atmospheric pressure	70КРа ~ 106КРа		

## 8. Symbols and description

3	Refer to instruction manual/booklet		Manufacturer
LOT	Batch code	$\sim$	Date of manufacture
SN	Serial number	$\square$	Use-by date
Ŕ	Type BF applied part	$\triangle$	Caution
*	Keep away from sunlight	ľ	Temperature limitation
) E	Humidity limitation	NON	Non-sterile
(1 <b>1</b> )	Single patient multiple use	Ť	Keep away from rain
<b>C €</b> <sub>0123</sub>	This device complies with the requirements of Regulation (EU) 745/2017	QTY	Quantity
<u>tt</u>	This way up	EC REP	Authorized representative
Ţ	Fragile; handle with care	(† (†	Atmospheric pressure limitation
X	The symbol indicating separate collection for EEE	MD	Medical Device

## 9. Contact information

Vega Technologies Inc.

No.505, Shengping Mid.Road, Yangwu, Dalang Town, 523789 Dongguan City, Guangdong Province, PEOPLE'S REPUBLIC OF CHINA



Website: www.vegamedical.com



MedPath GmbH Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany



Made by: 一刻 1 步

Approved by: 法家庭